


HIRANUMA APPLICATION DATA	Automatic Titrator	Data No.	A10	Nov 27, 2023
FOOD	Peroxide value measurement for cooking oil			

1. Abstract

Degraded fats and oils include peroxide chemical species. Therefore, determination of peroxide value is an effective analytical method to know the index of fat and oil degradation. This method is described in a variety of official standards in fats, oils, health science and pharmaceutical area.

Definition of peroxide value depends on respective standard method. For example, “Method of Analysis in Health Science” in Japan defines peroxide value as “milliequivalents of iodine generated from potassium iodide in 1 kg of sample [meq/kg]” (Formula (1)). The generated iodine is finally titrated with sodium thiosulfate (Formula (2)).

Example of titration for peroxide value in cooking oil is introduced here. The experimental procedure conforms to the “Japanese Agricultural Standards”.



2. Configuration of instruments and Reagents

(1) Instruments

- Main unit : Automatic Titrator COM series
- Electrode : Platinum-Reference combination electrode(long type) PR-733BZ
- Options : Shaft modifications to the stirrer supplied with the unit, buret tip modifications and modifications to the tube holder and other ancillary components are required. (Refer to “5. Note (4)”)

(2) Regents

- Titrant : 0.01 mol/L Sodium thiosulfate standard solution
- Additive : Saturated potassium iodide solution
Dissolve 80 g of potassium iodide in 50 mL of DI water.
- Solvent : Mixed solvent of 2,2,4-trimethylpentane and acetic acid with 2 : 3 ratio [v/v]
- Purge gas : Nitrogen gas

3. Measurement procedure

- (1) Take about 10 g of the sample into 200 mL Erlenmeyer flask with stopper and weigh it accurately.
- (2) Add 60 mL of mixed solvent to dissolve sample.
- (3) The air in flask is purged with nitrogen gas for 1 minutes. The flow rate of nitrogen gas is 300 mL/min.
- (4) Add 1 mL of saturated potassium iodide solution with pipette. Close the flask with stopper and swirl gently for 1 minute.
- (5) Leave the flask to stand for 5 minutes in dark place.
- (6) Add 60 mL of DI water to the flask.
- (7) Immerse the electrode and titrate with 0.01 mol/L sodium thiosulfate standard solution.
- (8) Measure the blank value by testing of (2) ~ (7) without sample.
- (9) After the measurement finished, the electrode should be immersed in mixed solvent, and additionally washed with DI water. Besides immerse it in DI water between each measurement to activate electrode.

4. Measurement conditions and Results

Examples of titration conditions

Measurement of blank


Condition No. 1			
Method	Auto	Constant No.	1
Buret No.	1	Size	0.00000 g
Amp No.	1	Blank	0.0000 mL
D.Unit	mV	Molality	0.0100 mol/L
S-Timer	30 sec	Factor	1.00400
CP mL	0 mL	K	0.00000
T-Timer	0 sec	L	0.00000
D.P. mL	0 mL	Unit	mL
End Sens	200	Formula	D
Over mL	0 mL	Digits	3
Max. Vol.	20 mL	Auto input Parameter	None
		Mode No.	15
		Pre Int	1 sec
		Del K	3
		Del Sens	0 mV
		Int time	5 sec
		Int Sens	3 mV
		Buret Speed	2
		Pulse	40
			0.05 mL

Measurement of sample


Condition No. 2			
Method	Auto	Constant No.	2
Buret No.	1	Size	0.00000 g
Amp No.	1	Blank	0.1270 mL
D.Unit	mV	Molality	0.0100 mol/L
S-Timer	30 sec	Factor	1.00400
CP mL	0 mL	K	10.00000
T-Timer	0 sec	L	0.00000
D.P. mL	0 mL	Unit	meq/kg
End Sens	200	Formula	(D-B)*K*F/S
Over mL	0.2 mL	Digits	4
Max. Vol.	20 mL	Auto input Parameter	None
		Mode No.	20
		Pre Int	1 sec
		Del K	5
		Del Sens	0 mV
		Int time	5 sec
		Int Sens	3 mV
		Buret Speed	2
		Pulse	40
			0.050 mL

Measurement results

Sample	Measurement	Size (g)	Titrant volume (mL)	Peroxide value (meq/kg)	Statistical result		
	No.						
BLANK	1	-	0.122	-	Avg.	0.127	mL
	2	-	0.131	-			
SAMPLE	1	10.0262	0.954	0.8281	Avg.	0.82	meq/kg
	2	10.0213	0.929	0.8035	SD	0.01	meq/kg
	3	10.0051	0.949	0.8219	RSD	1.56	%



Measurement of blank



Measurement of sample

Examples of measurement curves

5. Note

- (1) Please use Erlenmeyer flask with stopper to avoid sublimation of iodine and entraining of oxygen from air.
- (2) Timely measure the blank because saturated potassium iodide solution is easy to degenerate. Additionally, please periodically prepare fresh solution.
- (3) The type of solvent used in this measurement, the amount of solvent added, and the time left in the dark vary with each standard test methods. Please make the appropriate modifications according to the standard test methods used.
- (4) In this report, since a 200 mL Erlenmeyer flask is used as a titration vessel, platinum-reference combination electrode PR-733BZ (Long type) was used. The stirrer is modified to accommodate Erlenmeyer flasks as shown in Fig. 1 to attach the long type electrode.




Figure 1 A set of electrode and stirrer parts for Erlenmeyer flasks

Keywords : Cooking oil, Peroxide value, Redox titration