

HIRANUMA APPLICATION DATA	Automatic Titrator	Data No.	H4	Apr. 13, 2022
SODA PULP INDUSTRY	Determination of available chlorine in sodium hypochlorite			

1. Abstract

Sodium hypochlorite is used for bleaching and sterilization of tap water because of its strong oxidizing and disinfecting properties. Sodium hypochlorite is relatively stable at the alkaline region. However, it is unstable at the acidic region and becomes hypochlorous acid (HClO). It oxidizes water and generates chlorine (Cl₂). The concentration of available chlorine has to be measured regularly because sodium hypochlorite degrades slowly and generates sodium chloride.

This report introduces an example for determination of available chlorine as follows:

- 1) Add potassium iodide to sodium hypochlorite to generate free iodine.
- 2) Titrate the free iodine generated from the reaction (1) with sodium thiosulfate to determine available chlorine (2) by redox titration.



2. Configuration of instruments and reagents

(1) Configuration

Main unit : Hiranuma Automatic Titrator COM series
 Electrodes : Platinum electrode PT-301
 : Reference electrode RE-201Z

※Instead of above electrodes, the following electrodes are usable.

- PR-701BZ (Platinum reference combination electrode)
- Combination of PT-301 (Platinum Electrode) and GR-501BZ (Glass reference electrode)

<Remark> Measurement of alkaline component in sodium hypochlorite will be possible with the combination of PT-301 and GR-501BZ.

(2) Regents

Titrant : 0.3 mol/L Sodium thiosulfate standard solution
 Additive solution : 10 mL of 20% potassium iodide solution
 Buffer solution : 10 mL of 50 % acetic acid solution


3. Measurement procedure

- (1) Dispense 2 mL of sample and accurately weigh it.
- (2) Add about 10 mL of 20 % potassium iodide solution.
- (3) Add about 30 mL of DI water.
- (4) Add about 10 mL of 50 % acetic acid solution.
- (5) Immerse electrodes and titrate with 0.3 mol/L sodium thiosulfate standard solution.

4. Measurement conditions and results

Example of titration conditions

Cnd. No.	1	Constant No.	1	Mode No.	5
Method	Auto	Size	0 g	Pre Int	0 sec
Buret No.	1	Blank	0 mL	Del K	5
Amp No.	2	Morality	0.3 mol/L	Del Sens	0 mV
D.Unit	mV	Factor	1.013	Int Time	3 sec
S- Timer	15 sec	K	35.45	Int Sens	3 mV
C.P. mL	0 mL	L	0	Brt Speed	2
T.Timer	0 sec	Unit	%	Pulse	40
D.P. mL	0 mL	Formula	$(D-B)*K*F*M/(S*10)$		0.05 mL
End Sens	200	Decimal Places	4		
Over mL	1 mL	Auto input parameter	None		
Max Vol.	20 mL				



Example of titration curve

Measurement results

Number of Measurement	Size (g)	Titrant Volume (mL)	Concentration (%)
1	2.6123	30.710	12.6650
2	2.5984	30.528	12.6572
3	2.5988	30.572	12.6735
Statistic calculation	Avg.		12.67 %
	SD		0.01 %
	RSD		0.06 %

5. Note

Please refer to the following points to improve the measurement accuracy.

- 1) Weigh a sample size quickly because available chlorine is unstable.
- 2) Titrate free iodine immediately because the iodine generated by addition of potassium iodide readily volatilizes.

Addition of plenty potassium iodide is required because the volatilization of iodine depends on the concentration of the added potassium iodide.

Keywords: Sodium hypochlorite, Available chlorine, Redox titration

*Some measurement would not be possible depending on optional configuration of system.