



Inorganic acids &  
Mixed acidsFractional determination of  
hydrochloric acid and sulfuric acid

## 1. Abstract

The mixed solution of hydrochloric acid and sulfuric acid works as strong acid, and also has the strong oxidizability and solvency. It is used as the surface treatment solution for metals, glass products, and semiconductors. Hydrochloric acid and sulfuric acid are strong acids, therefore the fractional determination by neutralization titration is difficult.

The total acids of hydrochloric acid and sulfuric acid in the mixed solution are determined first by neutralization titration in this report. Nitric acid is added to the sample solution continuously to adjust the pH. Finally, the concentration of hydrochloric acid is determined by precipitation titration, the concentration of sulfuric acid is calculated by subtracting the concentration of hydrochloric from the total acids concentration.

The example of fractional determination for hydrochloric acid and sulfuric acid with additional burets are introduced in this report.



## 2. Configuration of instruments and reagents

## (1) Configuration of instruments

Main unit	:	Hiranuma Automatic Titrator COM series
Options	:	One buret, Peristaltic pump type dispenser
Electrodes	:	Glass electrode GE-101B for total acids measurement Connect to IE-1. Silver -reference electrode AGR-811Z for hydrochloric acid measurement (Double junction type) Connect to IE-2 and RE-2.

## (2) Reagents

- Titrant : 0.1 mol/L Sodium hydroxide standard solution for total acids measurement  
0.1 mol/L Silver nitrate for hydrochloric acid measurement
- Additive solution : 2 mL of 1 mol/L Nitric acid solution for pH adjustment

## 3. Measurement procedure

- (1) Dispense 1 mL of sample into 30 mL DI water in a volumetric flask with volumetric pipette. Cool it to room temperature and dilute to 100 mL with DI water.
- (2) Dispense 10 mL of diluted sample solution into a 100 mL beaker with volumetric pipette.
- (3) Add 40 mL of DI water.
- (4) Immerse electrodes and start titration with 0.1 mol/L sodium hydroxide standard solution.
- (5) After the procedure (4), 2 mL of 1 mol/L nitric acid is automatically dispensed by optional dispenser.
- (6) Titrate with 0.1 mol/L silver nitrate standard solution with optional buret.

## 4. Measurement conditions and results

### Examples of titration conditions

#### (1) Titration for total acids with sodium hydroxide standard solution (converted into hydrochloric acid)

Cndt No.	1	ConstantNo.	1	Mode No.	20
Method	Auto	Size	10.000 mL	Pre Int	0 sec
Buret No.	1	Blank	0 mL	Del K	9
Amp No.	1	Molarity	0.100 mol/L	Del Sens	0 mV
D. Unit	pH	Factor	1.004	Int Time	2 sec
S-Timer	5 sec	K	1.000	Int Sens	3 mV
C.P. mL	0 mL	L	0.010	Brst Speed	2
T Timer	0 sec	Unit	mol/L	Pulse	40
D.P. mL	0 mL	Formula	(D-B)*K*F*M/(S*L)		
End Sens	1000	Decimal Places	3		
Over mL	0 mL	Auto In Pram.	Non		
Max.Vol.	20 mL				

K: Equivalent of NaOH to HCl L: Dilution ratio

#### (2) Dispense 1 mol/L nitric acid

Cndt No.	2
Method	Disp
Buret No.	3
S-Timer	0 sec
Disp Vol.	2 mL

#### (3) Titration for hydrochloric acid with silver nitrate standard solution


Cndt No.	3	ConstantNo.	3	Mode No.	6
Method	Auto	Size	10.000 mL	Pre Int	0 sec
Buret No.	2	Blank	0.0000 mL	Del K	2
Amp No.	2	Molarity	0.100 mol/L	Del Sens	0 mV
D. Unit	mV	Factor	1.005	Int Time	3 sec
S-Timer	5 sec	K	1.000	Int Sens	3 mV
C.P. mL	0 mL	L	0.010	Brst Speed	2
T Timer	0 sec	Unit	mol/L	Pulse	40
D.P. mL	0 mL	Formula	(D-B)*K*F*M/(S*L)		
End Sens	300	Decimal Places	4		
Over mL	0.2 mL	Auto In Pram.	Non		
Max.Vol.	20 mL				

(4) Calculation for sulfuric acid concentration


Cndt No.	4	ConstantNo.	4
Method	Calc	Size	10.000 mL
		Blank	0.000 mL
		Molarity	0.100 mol/L
		Factor	1.000
		K	1.000
		L	0.01
		Unit	mol/L
		Formula	(CA-CB)/2
		Decimal Places	4
		Auto In Pram.	Non

Measurement results

Number of Measurement	Size (mL)	Dilution ratio	NaOH Titrant Volume (mL)	Total acids (mol/L) (converted into hydrochloric acid)	AgNO <sub>3</sub> Titrant Volume (mL)	Hydrochloric acid Conc. (mol/L)	Sulfuric acid Conc. (mol/L)
1			6.777	6.8041	0.683	0.6864	3.0589
2	10	1/100	6.776	6.8031	0.680	0.6834	3.0599
3			6.779	6.8061	0.682	0.6854	3.0603
Statistic calculation	Avg.			6.804 mol/L		0.685 mol/L	3.060 mol/L
	SD			0.002 mol/L		0.002 mol/L	0.001 mol/L
	RSD			0.022 %		0.223 %	0.02 %



Measurement of total acids



Measurement of hydrochloric acid

Examples of titration curves

## 5. Note

### (1) Fractional determination of hydrochloric acid and sulfuric acid

The concentration of hydrochloric acid is measured by the determination of chloride ion in this method.

This method cannot be applied to the sample containing chloride ion from other than hydrochloric acid.

### (2) Control of titrant

The carbon dioxide gas absorber (soda lime) on reagent bottle has to be regularly exchanged because sodium hydroxide used for total acids determination readily absorbs carbon dioxide gas in the air.

Keywords: Fractional determination of hydrochloric acid and sulfuric acid, Neutralization titration,  
Precipitation titration