

HIRANUMA APPLICATION DATA	Automatic Titrator	Data No.	D2	Nov. 14, 2018
Environment	Total hardness in tap water			

1. Abstract


Hardness in water is total amount of calcium and magnesium ion in water converted to mg/L of comparable calcium carbonate. “Hardness” has the following types:

- (1) Total hardness Total amount of calcium and magnesium ions
- (2) Calcium hardness Calcium ion
- (3) Magnesium hardness (Total hardness) – (Calcium hardness)
- (4) Non-carbonate hardness (Permanent hardness)
- (5) Carbonate hardness (Temporary hardness)

These are stipulated in some applicable standard such as *Standard Methods for the Examination of Water*, *Standard Methods of Analysis for Hygienic Chemists*, and *JIS K0101 Testing methods for industrial water* etc.

This report introduces an example for determination of total hardness in tap water with photometric titration method using EDTA standard solution according to the *Standard Methods for the Examination of Water*.

Take 100 mL of sample and add 1 mL of 0.01 mol/L magnesium chloride, 2 mL ammonia buffer, and 0.2 mL of EBT indicator. Titrate with 0.01 mol/L EDTA standard solution (red → blue color). Perform the same procedure for 100 mL of DI water instead of sample as blank measurement.



2. Configuration of instruments and reagents

(1) Configuration of instruments

Main unit : Hiranuma Automatic Titrator COM Series
(M type photometric unit for photometric titration with 650 nm optical filter)

(2) Reagents

Titrant : 0.01 mol/L EDTA standard solution
Additive solution : 0.01 ml/L Magnesium chloride solution
Buffer solution : Ammonium chloride-ammonia buffer solution (pH 10)
Dissolve 6.75 g of ammonium chloride in 30 mL of DI water, add 57 mL of ammonia solution, and dilute to 100 mL with DI water.
Indicator : EBT Indicator
Dissolve 0.5 g of eriochrome black T and 4.5 g of hydroxylammonium chloride in 100 mL of ethanol (95 % [v/v]).

3. Measurement procedure

- (1) Dispense 100 mL of sample into a 200 mL beaker with volumetric pipette.
- (2) Add 1 mL of 0.01 mol/L magnesium chloride solution with volumetric pipette.
- (3) Add 2 mL of ammonia buffer.
- (4) Add 0.2 mL of EBT indicator.
- (5) Immerse photometric probe and start titration with 0.01 mol/L EDTA standard solution.
- (6) Perform the same procedure for 100 mL of DI water as blank measurement

4. Measurement conditions and results



Example of titration condition

Measurement of blank

Cndt No.	1	ConstantNo.	1	Mode No.	8
Method	B cross	Size	0 mL	Pre Int	0 sec
Buret No.	1	Blank	0 mL	Del K	5
Amp No.	2	Molarity	0.01 mol/L	Del Sens	0 mV
D. Unit	T%	Factor	1.001	Int Time	5 sec
S-Timer	10 sec	K	0	Int Sens	3 mV
C.P. mL	0 mL	L	0	Brt Speed	2
T Timer	0 sec	Unit	mL	Pulse	40
D.P. mL	0 mL	Formula	D		
End Sens	500	Digits	4		
Over mL	1 mL	Auto In Pram.	Non		
Max.Vol.	20 mL				

Measurement of sample

Cndt No.	2	ConstantNo.	2	Mode No.	8
Method	B cross	Size	100 mL	Pre Int	0 sec
Buret No.	1	Blank	1.040 mL	Del K	5
Amp No.	2	Molarity	0.01 mol/L	Del Sens	0 mV
D. Unit	T%	Factor	1.001	Int Time	5 sec
S-Timer	10 sec	K	100	Int Sens	3 mV
C.P. mL	0 mL	L	0	Brt Speed	2
T Timer	0 sec	Unit	PPM	Pulse	40
D.P. mL	0 mL	Formula	$(D-B)*K*F*M*1000/S$		
End Sens	500	Digits	4		
Over mL	1 mL	Auto In Pram.	Non		
Max.Vol.	20 mL				



Examples of titration curves

Measurement results

Measurement of blank

Number of Measurement	Size (mL)	Titrant Volume (mL)
1		1.041
2	-	1.041
3		1.038
Statistic calculation	Avg.	1.040 mL
	SD	0.0017 mL
	RSD	0.17 %

Measurement of sample

Number of Measurement	Size (mL)	Titrant Volume (mL)	Total hardness (ppm)
1	100	6.718	56.837
2	100	6.745	57.107
3	100	6.719	56.847
Statistic calculation		Avg.	56.9 ppm
		SD	0.153 ppm
		RSD	0.27 %

5. Note

(1) Measurement

The total hardness in tap water can accurately be determined by photometric titration method with EBT indicator using photometric probe.

(2) End point detection

“Method” on the condition parameter is set to “B cross” because the point where indicator color change is completed is detected as endpoint. When the concentration of total hardness is relatively high and much titrant volume is required, “C.P. mL” function is useful. Titrant can be added at once by setting “C.P. mL” to a few mL smaller than titrant volume at the endpoint, it can reduce measurement time.

Keywords: Tap water, Total hardness, Photometric titration, EDTA standard solution